EASYspin Universal Plant RNA Kit**EASYspin 通用植物 RNA 快速提取试剂盒/DNase I**

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1. 番茄叶: Effect of Low Temperature Stress on the Expression of ProDH Gene and the Activities of the Proline Dehydrogenase in Leaves of Tomato Seedling. *Chinese Agricultural Science Bulletin* 2012,28(10):132-135
2. 大豆: RNA-seq Analysis Reveals Ethylene-Mediated Reproductive Organ Development and Abscission in Soybean(*Glycine max* L. Merr.). *Plant Mol Biol Rep*, 2012, published online: 4 Dec, 2012
3. 大豆 1: Construction of ethylene regulatory network based on the phytohormones related gene transcriptome profiling and prediction of transcription factor activities in soybean. *Acta Physiol Plant*, 2012, published online: 12 Dec, 2012
4. 红花玉兰: Expression Analysis of MAwuAG in Different Organs and Developmental Stages of *Magnolia wufengensis*. *Chinese Bulletin of Botany*, 2013, 48 (2): 1-5
5. 毛桃: Cloning and Phylogeny Analysis of PpAP2 Floral Homologous Genes in Peach. *Chinese Agricultural Science Bulletin*, 2013, 29(7): 99-104
6. 茶梅花瓣: Comparison and Analysis of Methods of Extracting Total RNA from Petals of *Camellia sasanqua*. *Chinese Agricultural Science Bulletin*.2013,29(28):129-133.
7. 棉花: Analysis of sea-island cotton and upland cotton in response to *Verticillium dahliae* infection by RNA sequencing. Sun et al. *BMC Genomics* 2013, 14:852 /1471-2164/14/852.
8. 荞麦和拟南芥: Ectopic expression of *FaesAP3*, a *Fagopyrum esculentum* (Polygonaceae) AP3 orthologous gene rescues stamen development in an *Arabidopsis ap3* mutant. *Gene* 2014, 550(2): 200-206
9. 玫瑰花: Precise spatio-temporal modulation of ACC synthase by MPK6 cascade mediates the response of rose flowers to rehydration. *The Plant Journal* 2014, 79(6): 941-950
10. 棉花和拟南芥: Functional characterization of GhAKT1, a novel Shaker-like K⁺ channel gene involved in K⁺ uptake from cotton (*Gossypium hirsutum*). *Gene* 2014, 545(1): 61-71
11. 棉花 1: Gibberellin Overproduction Promotes Sucrose Synthase Expression and Secondary Cell Wall Deposition in Cotton Fibers. *PLoS ONE* 2014, 9(5): e96537. doi:10.1371/journal.pone.0096537
12. 苹果: Low Medium pH Value Enhances Anthocyanin Accumulation in *Malus Crabapple* Leaves. *PLoS ONE* 2014, 9(6): e97904. doi:10.1371/journal.pone.0097904